

Neurotone AI

# Integrating Sound Therapy and Cognitive Behavioral Intervention in Tinnitus Care



Rick Carlson, CEO




Dr. Miles Aron, COO

# What's new with Lace

- 1k+ practice partners in 18 months
- Dr. Brian Taylor – VP of Clinical Research and Professional Relations
  - Editor of *Audiology Practices*, quarterly ADA journal
- Coming soon: Targeted Aural Rehab
- Now available:
  - Train in 10 languages
  - Provider Prescribed Programs

Je m'appelle Lacey,  
je suis votre coach  
IA. C'est parti!

  
Bienvenu!





# We're built on a foundation of science




**Jennifer Henderson Sabes, AuD**  
Co-inventor & Scientific Advisor

Senior Research Audiologist @ UCSF  
Audiology Professor @ Pacific  
AAA fellow and ASHA CCC-A certified



**Robert Sweetow, PhD**  
Co-Inventor & Scientific Advisor

Professor Emeritus of Otolaryngology @ UCSF  
Former Chief of Audiology @ UCSF  
AAA Distinguished Achievement Award



**Miles Aron, PhD**  
Co-Founder & COO

PhD in Biomedical Acoustics @ Oxford  
US Fulbright Scholar  
BSE in Acoustical Engineering @ Hartford



**Nancy Tye-Murray, PhD**  
Scientific Advisor

Professor of Otolaryngology  
@ Washington University  
School of Medicine



# The Science

# Habituation Therapy

Based on Tinnitus Retraining Therapy (TRT)

5 studies · 541 participants

KEY FINDING	STUDY	METHOD	N
<b>+54% ↑</b> THI improvement at 6 mo.	<b>Patient-Based Outcomes In Patients With Primary Tinnitus Undergoing TRT</b> <i>Berry et al., JAMA Otolaryngology, 2002</i>	Prospective, nonrandomized; THI before & after TRT	n=32
<b>82% ↓</b> Showed 20-point THI decrease	<b>Implementing The Neurophysiological Model Of Tinnitus Retraining Therapy</b> <i>Jastreboff &amp; Hazell, 2004</i>	Counseling + sound therapy based on neurophysiological model	n=303
<b>72% ↓</b> Significant decrease at 12 mo., Improvements over time support multiple practical delivery pathways	<b>Open Ear Hearing Aids In Tinnitus Therapy Vs. Sound Generators</b> <i>Parazzini et al., Int. J. Audiology, 2011</i>	Randomized prospective parallel-group study	n=91
<b>TFI ↑</b> 33 pts TFI improvement (HA+), 21 pts with HA, 31 pts with extended-wear HA	<b>Rct Comparing Extended-Wear, Conventional Has, And Combination Instruments</b> <i>Henry et al., Journal of AAA, 2016</i>	RCT comparing 3 bilateral hearing device types	n=55
<b>60 studies</b> Tinnitus distress consistently associated with certain traits	<b>Anxiety, Depression, And Personality Traits Relevant To Tinnitus</b> <i>Durai &amp; Searchfield, Frontiers of Neuroscience, 2017</i>	Scoping review of the literature	n=60

# Adaptation Therapy

4 studies · 174 participants

Notched audio designed to drive cortical rebalancing (lateral inhibition + plasticity)

KEY FINDING	STUDY	METHOD	N
<h2>Tinnitus ↓↓</h2> <p>Tinnitus loudness and tinnitus-related auditory cortex activity after 12 mo.</p>	<p>Tailor-Made Notched Music Reduces Tinnitus Loudness and Auditory Cortex Activity</p> <p><i>Okamoto et al., PNAS, 2010</i></p>	<p>12-mo double-blind longitudinal study with target, placebo, and monitoring groups</p>	<p>n=23</p>
<h2>At 3 months</h2> <p>TMNMT appeared more effective than TRT</p>	<p>Efficacy of Tailor-Made Notched Music Training Vs. TRT in Adults With Chronic Subjective Tinnitus</p> <p><i>Tong et al., Ear &amp; Hearing, 2023</i></p>	<p>Single-blind randomized controlled trial</p>	<p>n=112</p>
<h2>Outcomes vary</h2> <p>By stimulus and patient factors</p>	<p>Comparison of Acoustic Therapies for Tinnitus Suppression: A Preliminary Trial</p> <p><i>Schad et al., Int. J. Audiology, 2017</i></p>	<p>Evaluated notched noise against other sound controls</p>	<p>n=30</p>
<h2>Tinnitus ↓↓</h2> <p>Reduced tinnitus loudness and tinnitus-related neural activity after 3 days, effects accumulated across sessions</p>	<p>Inhibition-Induced Plasticity in Tinnitus Patients After Repetitive Exposure to Notched Music</p> <p><i>Pantev et al., BMC Neurology, 2016</i></p>	<p>3-month protocol in chronic tonal tinnitus using tailor-made notched music (2 hrs/day)</p>	<p>n=9</p>

# Relief Mode

Short partial–masking sessions for spikes + residual inhibition potential

3 studies · 361 participants

KEY FINDING	STUDY	METHOD	N
<b>TFI ↑↑</b> Significant improvement in tinnitus–related distress (TFI), supported residual inhibition technique feasibility	<b>Evaluating a Residual Inhibition Technique in Hearing Aids for Tinnitus Suppression</b> <i>Quinn et al., Am J Audiol, 2023</i>	Pilot hearing aid sound therapy with 2–month intervention and 1–month washout	n=20
<b>20 points ↓</b> THI Improvement With Amplification + Masking	<b>Tinnitus Sound Amplification and Masking Therapy: Outcomes at UNC Chapel Hill</b> <i>Meade et al., Am J Otolaryngology, 2024</i>	Review of tinnitus patients treated over 12 months	n=141
<b>~50%</b> Experienced Residual Inhibition when using Masking	<b>Management of Tinnitus: Discussion Paper</b> <i>Hazell, Journal of the Royal Academy, 1985</i>	Discussion paper based on tinnitus clinic experience	n=200

# Sleep Mode

Maximizes sound therapy time when patients' minds are most open to habituation

2 studies · 69 participants

KEY FINDING	STUDY	METHOD	N
<h2>Tinnitus ↓↓</h2> <p>Significant reduction in intensity &amp; burden after 3 months</p>	<h3>Impact of Reduction of Tinnitus Intensity on Patients' Quality of Life</h3> <p><i>Drexler et al., Int. J. Audiology, 2015</i></p>	Study of nightly customized sound stimulation during sleep	n=11
<h2>62% ↓</h2> <p>Average reduction in perceived tinnitus and 77% average improvement in TFI when using Personalized sound therapy vs. generic sounds</p>	<h3>RCT of a Novel Device for Tinnitus Sound Therapy During Sleep</h3> <p><i>Theodoroff et al., Am J Audiology, 2017</i></p>	RCT of 3 sleep-based sound therapy approaches	n=58

# CBT, Meditation and Mindfulness Tools

3 studies · 170 participants

Helps patients identify and restructure unhelpful thoughts, build coping strategies, and improve function


KEY FINDING	STUDY	METHOD	N
<h2>MBCT wins</h2> <p>MBCT outperformed relaxation training in reducing tinnitus severity, with lasting benefit at 6 months</p>	<p>Mindfulness-Based Cognitive Therapy for Chronic Tinnitus: A Randomized Controlled Trial</p> <p><i>McKenna, Marks &amp; Schaette, Psychother Psychosom, 2017</i></p>	<p>RCT comparing 8 weekly 120-min MBCT sessions with intensive relaxation training</p>	<p>n=75</p>
<h2>Meditation wins</h2> <p>Mindfulness meditation produced greater improvement than relaxation therapy, and a larger reduction in tinnitus distress</p>	<p>A Randomised Controlled Study of Mindfulness Meditation Vs. Relaxation Therapy in the Management of Tinnitus</p> <p><i>Arif et al., Laryngol Otol, 2017</i></p>	<p>RCT of 5-session mindfulness vs. relaxation treatment</p>	<p>n=86</p>
<h2>Wellbeing ↑↑</h2> <p>MBCT shifted resistance to acceptance</p>	<p>I Wasn't at War With the Noise: How MBCT Changes Patients' Experiences of Tinnitus</p> <p><i>Marks et al., Front Psychol, 2020</i></p>	<p>Qualitative follow-up within RCT; 5-session MBCT vs. relaxation</p>	<p>n=9</p>

# Introducing Tinnitus Pro

# Introducing Tinnitus Pro

The common sense tinnitus solution for the 90% of people that want a serious, science-based treatment at an affordable price.

- Best-in-class tinnitus treatment
- Gamified and engaging
- Customized for your practice
- Pro Portal patient management
- Personalized tinnitus tone finder
- Lifetime patient licenses



## Personalized Sound Therapy



## Cognitive Behavioral Intervention



CBI



Meditation





Available in 10+ languages

# Personalized Sound Therapy

# Therapeutic tone finder

First of its kind innovation in sound therapy

- Helps patients pinpoint their tinnitus tone with more precision for maximum benefit
- Personalized sound therapy is generated from the patient's unique tinnitus tone
- Patients experience pure tones, harmonics, and band noises all of which are captured here  
*(Tinnitus: Theory and Management, Henry, 2004.)*
- Powering the largest tinnitus study ever, advancing how we understand and treat tinnitus




Blend therapeutic tones  
with soundscapes

# Personalized sound therapy

Advanced care based on patients' own tinnitus tone

The only tinnitus app to provide patients with the option to treat tinnitus using therapies developed on the theories of neural habituation, adaptation, and inhibition.



## Why we've gone further

- Research suggests that different patients respond to different therapy modes
- Hearing care providers and their patients can find what works best




# Habituation Therapy

Therapeutic sounds with blendable soundscapes  
below masking level leads to tinnitus fading

Daily cumulative goal:

**1.5+ hours**   
More time can produce bigger gains



## Neural habituation

- Reduction in perceived loudness
- 25 pt reduction (54% improvement) in THI scores
- Key ingredient in Tinnitus Retraining Therapy

Berry et al, JAMA Otolaryngology, 2002.  
Jastreboff and Hazell, 2004.  
Parazzini et al. International Journal of Audiology, 2011.  
Henry et al. Journal of AAA 2016.  
Durai & Searchfield, Frontiers of Neuroscience, 2017.




# Adaptation Therapy

Adaptation therapy (aka “notched noise therapy”) removes your tinnitus frequency from therapeutic audio.

Daily goal:

**1.5+ hours** 

More time can produce bigger gains



Why it works?

## Neural adaptation

- Reduces hyperactivity through lateral inhibition and cortical plasticity (sound-excited neurons inhibit their overactive neighbors)
- 3-month RCT found greater improvements with adaptation vs. traditional sound therapy

Okamoto et al., PNAS 2010.

Tong et al., Ear & Hearing 2023.

Schad et al., International Journal of Audiology 2017.

Pantev et al., BMC Neurology 2016.

Confidential 17



# Relief Mode


Short partial-masking sessions with sounds based on the patient's unique tinnitus tone for instant relief

Daily goal:

**2-15 min**



During tinnitus spiking



Why it works?

- Puts patients in control, reduces contrast, attentional, and emotional factors
- 20 pt reduction in THI (52% improvement) in 2024 study
- ~50% of patients experience additional residual inhibition—a temporary post-sound quieting

Quinn et al., Am J Audiol 2023.

Meade et al., American Journal of Otolaryngology, 2024.

Hazell, Journal of the Royal Academy 1985.




# Sleep Mode

Maximizes time on sound therapy when patients' minds are most open to habituation

Daily goal:

**5–8 hours** 

Start at bedtime



## Why it works?

- Reduces tinnitus–silence contrast to ease sleep onset
- Personalized sound therapy outperforms generic sounds
- After 3 months: 62% average reduction in perceived tinnitus
- 77% average improvement in Tinnitus Functional Index (TFI)

Drexler et al, 2015 International Journal of Audiology  
Theodoroff et al, 2017 American Journal of Audiology

The background features a series of concentric, slightly irregular circles in a lighter shade of blue, creating a ripple effect that frames the central text.


# Cognitive Behavioral Intervention



# Cognitive behavioral intervention

50+ lessons from the leading audiology researchers behind LACE

Learn to identify and restructure negative thoughts to better cope with tinnitus.



Understanding Habituation

Reframing Negative Thoughts

Celebrate Small Wins

Use Sound During Quiet Tasks

Practice Self-Compassion

Review Medications

Reinforce Positive Moments

Try Sound Therapy Before Bed

Try a Mix of Tones and Noise

Stick to a Sleep Schedule

Use Sound Therapy at Night

Don't Lie Awake for Too Long

Try Short, Frequent Breaks


Be Curious, Not Critical

Try Deep Breathing



# Meditation and mindfulness tools

Guided sessions layered with personalized therapeutic sounds to reduce reactivity to tinnitus.



- Reduces tinnitus severity and distress
- Increases acceptance
- Reduces catastrophizing
- Reorients attention to neutral anchors
- Provides benefits long after meditating

Psychother Psychosom. 2017 (McKenna/Marks/Schaette)  
Laryngol Otol. 2017 (Arif et al.)  
Front Psychol. 2020 (Marks et al.)

How we help  
patients succeed





# Built for patient engagement

So that more patients stick with it and get results that last

Streaks help keep patients training their brain to tune out the tinnitus.

Avatar-based coaches keep patients engaged with CBI lessons.



Progress tracking so patients can stay informed and motivated to make progress.


Automated encouragement nudges patients toward therapy when they need it most and dials back as the tinnitus fades.



# Integrated surveys

Measure patient outcomes with more fidelity

- Tinnitus Handicap Inventory (THI)
- Tinnitus Functional Index (TFI)



Identify, quantify, and evaluate the difficulties that patients may be experiencing, and how well they are functioning in their daily lives because of tinnitus



# Pro-Portal: For professionals

## Patient management

- Automated patient messaging
- New patient sign up

**Tinnitus sessions**

Last session: Sept. 15, 2025  
Total session time: 8 hr, 23 min  
Longest session: 1 hr, 54 mins

**Total session time by mode**

Mode	Hours
Therapy	4:33:40
Relief	3:50:15
Sleep	0
Meditation	0

**Tinnitus perception**

**My progress**  
Track your tinnitus sessions and symptoms over time.

**Tinnitus sessions**  
Total hours: 59 hrs, 0 min

- 41:00:00 Tone targeting therapy time
- 41:00:00 Tone masking therapy time
- 4:00:00 Relief time

**Message from Dr. Malbaonta**

**Dr. Laura Jarman**  
1234 Main St  
Gainesville, FL

**19 days in a row!**  
Great job, keep going!

**Tinnitus check-in**  
Please answer 5 quick questions about how you've been feeling over the past week.

## Patient monitoring

- Patient progress
- Patient achievements
- Survey results (THI, TFI, Check-ins)



## Practice personalization

- Custom branding
- Hours and location
- Appointment scheduling



# Lace AI Pro + Tinnitus Pro

Better together

Patient gets a unified experience



Providers can prescribe one app or both



Provider gets all patient data in Pro Portal


# Lifetime licenses

Provide lifelong care for your patients

Tinnitus patients will unfortunately experience spikes over the course of months or years.

Why?

- Medications
- Lifestyle changes
- Deterioration, etc.



Our price

**TO PATIENTS MSRP \$499**

Partners set the price

All plans come with unlimited licenses.

# Patient Outcomes: First Look

# Early outcomes for bothersome tinnitus with Tinnitus Pro

April 2026 | 5 months post-launch | 25 patient subgroup focus

**40% ↑**

Clinically meaningful gains  
(>7pt THI reduction; Zeman 2011 )

**13 pts ↓**

Median THI change among  
improved patients

**19h**



Avg sound therapy  
per patient

SAMPLE SIZE

**100**





Patients who used the app and completed at least 2 surveys — the basis for this analysis

About 50% of patients with an incoming THI of 22+ improve

# Patient Spotlights

De-identified patients with clinically meaningful THI improvement

5 patients

PATIENT	THI CHANGE	PRIMARY MODE	FREQ
Patient A	 66 (Severe) → 32 (Mild) <b>-34 pts</b>	Habituation Therapy · 0 CBT videos	9,694 Hz
Patient B	 40 (Moderate) → 14 (Slight) <b>-26 pts</b>	Relief Mode · 23 CBT videos	4,677 Hz
Patient C	 50 (Moderate) → 30 (Mild) <b>-20 pts</b>	Relief Mode · 15 CBT videos	4,877 Hz
Patient D	 56 (Moderate) → 42 (Moderate) <b>-14 pts</b>	Habituation Therapy · 21 CBT videos	5,256 Hz




Auditory  
Nerve

Cochlea

# Let's recap

The common sense tinnitus solution for the 90% of people that want a serious, science-based treatment at an affordable price.

- Best-in-class tinnitus treatment
- Personalized tinnitus tone finder
- Gamified and engaging
- Customized for your practice
- Pro Portal patient management
- Lifetime patient licenses



## Personalized Sound Therapy



## Cognitive Behavioral Intervention



